

BLADEX, INC.

ADVANCED TECHNOLOGIES FOR MEDICAL DEVICES

**CHANGE YOUR BLADES
IN ONLY 3 SECONDS!**

INSERT HANDLE
WITH USED BLADE



TURN 90° COUNTER-CLOCKWISE
TO DISCHARGE THE USED BLADE



TO ATTACH A FRESH BLADE,
PUSH FORWARD AND TURN
90° CLOCKWISE



YOUR PRESCRIPTION
FOR SAFETY!



WWW.BLADEX.NET

